

52. van Noesel CJ, Lankester AC, van Schijndel GM, van Lier RA: The CR2/CD19 complex on human B cells contains the src-family kinase Lyn. *Int Immunol* 1993; 5: 699–705
53. Uckun FM, Burkhardt AL, Jarvis L, Jun X, Stealey B, Dibirdik I, Myers DE, Tuel-Ahlgren L, Bolen JB: Signal transduction through the CD19 receptor during discrete developmental stages of human B-cell ontogeny. *J Biol Chem* 1993; 268: 21172–21184
54. Chalupny NJ, Kanner SB, Schieven GL, Wee SF, Gilliland LK, Aruffo A, Ledbetter JA: Tyrosine phosphorylation of CD19 in pre-B and mature B cells. *Embo J* 1993; 12: 2691–2696
55. Bireland ML, Monroe JG: Biochemistry of antigen receptor signaling in mature and developing B lymphocytes. *Crit Rev Immunol* 1997; 17: 353–385
56. Fujimoto M, Fujimoto Y, Poe JC, Jansen PJ, Lowell CA, DeFranco AL, Tedder TF: CD19 is a specialized transmembrane adapter protein that amplifies Src-family protein tyrosine kinase activation. (submitted);
57. Buhl AM, Pleiman CM, Rickert RC, Cambier JC: Qualitative regulation of B cell antigen receptor signaling by CD19: Selective requirement for PI 3-kinase activation, inositol-1,4,5-trisphosphate production and Ca²⁺ mobilization. *J Exp Med* 1997; 186: 1897–1910
58. Li HL, Forman MS, Kurosaki T, Pure E: Syk is required for BCR-mediated activation of p90Rsk, but not p70S6k, via a mitogen-activated protein kinase-independent pathway in B cells. *J Biol Chem* 1997; 272: 18200–18208
59. El-Hillal O, Kurosaki T, Yamamura H, Kinet JP, Scharenberg AM: syk kinase activation by a src kinase-initiated activation loop phosphorylation chain reaction. *Proc Natl Acad Sci USA* 1997; 94: 1919–1924
60. Buhl AM, Cambier JC: Phosphorylation of CD19 Y484 and Y515, and linked activation of phosphatidylinositol 3-kinase, are required for B cell antigen receptor-mediated activation of Bruton's tyrosine kinase. *J Immunol* 1999; 162: 4438–4446
61. Tuveson DA, Carter RH, Soltoff SP, Fearon DT: CD19 of B cells as a surrogate kinase insert region to bind phosphatidylinositol 3-kinase. *Science* 1993; 260: 986–989
62. Weng WK, Jarvis L, LeBien TW: Signaling through CD19 activates Vav/mitogen-activated protein kinase pathway and induces formation of a CD19/Vav/phosphatidylinositol 3-kinase complex in human B cell precursors. *J Biol Chem* 1994; 269: 32514–32521
63. Li X, Sandoval D, Freeberg L, Carter RH: Role of CD19 tyrosine 391 in synergistic activation of B lymphocytes by coligation of CD19 and membrane Ig. *J Immunol* 1997; 158: 5649–5657
64. O'Rourke LM, Tooze R, Turner M, Sandoval DM, Carter RH, Tybulewicz VL, Fearon DT: CD19 as a membrane-anchored adaptor protein of B lymphocytes: costimulation of lipid and protein kinases by recruitment of Vav. *Immunity* 1998; 8: 635–645
65. Sato S, Jansen PJ, Tedder TF: CD19 and CD22 expression reciprocally regulates tyrosine phosphorylation of Vav protein during B lymphocyte signaling. *Proc Natl Acad Sci USA* 1997; 94: 13158–13162
66. Zhou LJ, Ord DC, Hughes AL, Tedder TF: Structure and domain organization of the CD19 antigen of human, mouse and guinea pig B lymphocytes. Conservation of the extensive cytoplasmic domain. *J Immunol* 1991; 147: 1424–1432
67. Zhou S, Cantley LC: Recognition and specificity in protein tyrosine kinase-mediated signalling. *Trends Biochem Sci* 1995; 20: 470–475
68. Clark MR, Johnson SA, Cambier JC: Analysis of Ig-alpha-tyrosine kinase interaction reveals two levels of binding specificity and tyrosine phosphorylated Ig-alpha stimulation of Fyn activity. *Embo J* 1994; 13: 1911–1919
69. Ruzzene M, Brunati AM, Donella-Deana A, Marin O, Pinna LA: Specific stimulation of c-Fgr kinase by tyrosine-phosphorylated (poly)peptides: possible implication in the sequential mode of protein phosphorylation. *Eur J Biochem* 1997; 245: 701–707
70. Tedder TF, Tuscano J, Sato S, Kehrl JH: CD22, a B lymphocyte-specific adhesion molecule that regulates antigen receptor signaling. *Annu Rev Immunol* 1997; 15: 481–504
71. Tooze RM, Doody GM, Fearon DT: Counterregulation by the coreceptors CD19 and CD22 of MAP kinase activation by membrane immunoglobulin. *Immunity* 1997; 7: 59–67
72. Weng WK, Shah N, O'Brien D, Van Ness B, LeBien TW: Differential induction of DNA-binding activities following CD19 cross-linking in human B lineage cells. *J Immunol* 1997; 159: 5502–5508
73. Li X, Carter RH: Convergence of CD19 and B cell antigen receptor signals at MEK1 in the ERK2 activation cascade. *J Immunol* 1998; 161: 5901–5908
74. Lankester AC, van Schijndel GM, Rood PM, Verhoeven AJ, van Lier RA: B cell antigen receptor cross-linking induced tyrosine phosphorylation and membrane translocation of a multimeric Shc complex that is augmented by CD19 co-ligation. *Eur J Immunol* 1994; 24: 2818–2825
75. Beckwith M, Jorgensen G, Longo DL: The protein product of the proto-oncogene c-cbl forms a complex with phosphatidylinositol 3-kinase p85 and CD19 in anti-IgM-stimulated human B-lymphoma cells. *Blood* 1996; 88: 3502–3507
76. Chalupny NJ, Aruffo A, Esselstyn JM, Chan PY, Bajorath J, Blake J, Gilliland LK, Ledbetter JA, Tepper MA: Specific binding of Fyn and phosphatidylinositol 3-kinase to the B cell surface glycoprotein CD19 through their src homology 2 domains. *Eur J Immunol* 1995; 25: 2978–2984
77. Marte BM, Downward J: PKB/Akt: connecting phosphoinositide 3-kinase to cell survival and beyond. *Trends Biochem Sci* 1997; 22: 355–358
78. Franke TF, Kaplan DR, Cantley LC: PI3K: downstream AKTion blocks apoptosis. *Cell* 1997; 88: 435–437
79. Suzuki H, Terauchi Y, Fujiwara M, Aizawa S, Yazaki Y, Kadowaki T, Koyasu S: Xid-like immunodeficiency in mice with disruption of the p85 α subunit of phosphoinositide 3-kinase. *Science* 1999; 283: 390–392
80. Fruman DA, Snapper SB, Yballe CM, Davidson L, Yu JY, Alt FW, Cantley LC: Impaired B cell development and proliferation in absence of phosphoinositide 3-kinase p85 α . *Science* 1999; 283: 393–397
81. Fujimoto M, Bradney AP, Poe JC, Steeber DA, Tedder TF: Modulation of B lymphocyte antigen receptor signal transduction by a CD19/CD22 regulatory loop. *Immunity* 1999; 11: 191–200
82. Engel P, Nojima Y, Rothstein D, Zhou LJ, Wilson GL, Kehrl JH, Tedder TF: The same epitope on CD22 of B lymphocytes mediates the adhesion of erythrocytes, T and B lymphocytes, neutrophils and monocytes. *J Immunol* 1993; 150: 4719–4732
83. Kelm S, Pelz A, Schauer R, Filbin MT, Tang S, de Bellard ME, Schnaar RL, Mahoney JA, Hartnell A, Bradfield P, Crocker PR: Sialoadhesin, myelin-associated glycoprotein and CD22 define a new family of sialic acid-dependent adhesion molecules of the immunoglobulin superfamily. *Curr Biol* 1994; 4: 965–972
84. Schulte RJ, Campbell MA, Fischer WH, Sefton BM: Tyrosine phosphorylation of CD22 during B cell activation. *Science* 1992; 258: 1001–1004
85. Wilson GL, Fox CH, Fauci AS, Kehrl JH: cDNA cloning of the B cell membrane protein CD22: a mediator of B-B cell interactions. *J Exp Med* 1991; 173: 137–146